

# Complete FLASH Installation Protocol

## Muscle Fiber Analyzer - Complete Installation

### System Requirements

#### Minimum configuration required:



- **Operating System:** Windows 10/11 (x64), macOS 10.15+, or Linux Ubuntu 18.04+
  - **RAM:** 8 GB minimum (16 GB recommended)
  - **Disk Space:** 5 GB free space
  - **Graphics Card:** NVIDIA (optional, for GPU acceleration)
  - **Administrator Rights:** Required only for GPU/CUDA installation
- 

## STEP 0: GPU Verification (Windows only - optional)

### A. Graphics Card Identification

1. **Open Task Manager:**
  - Press `Ctrl + Shift + Esc`
  - Click on the "Performance" tab
  - Select "GPU" in the left panel
2. **Check your graphics card:**
  - **Required:** NVIDIA GTX 10xx, RTX series, or newer
  - **Recommended GPU RAM:** 6 GB or more
  - Note the exact card name (e.g., "NVIDIA GeForce RTX 3070")

#### Compatibility check:

-  **Compatible:** GTX 1060, RTX 2060, RTX 3070, RTX 4080, etc.
-  **Not compatible:** GTX 900 series or older, AMD cards

### B. CUDA Toolkit 11.8 Installation (Administrator rights required)

1. **Download CUDA Toolkit 11.8:**
  - Go to: <https://developer.nvidia.com/cuda-11-8-0-download-archive>
  - Select: Windows → x86\_64 → 10/11 → exe (local)
2. **CUDA Installation:**
  - Close all applications
  - Run the downloaded file as administrator
  - Choose "Custom installation" if offered
  - Restart computer after installation

### 3. CUDA Verification:

- Open command prompt (cmd)
  - Type: `nvcc --version`
  - Verify that "Cuda compilation tools, release 11.8" is displayed
- 

## STEP 1: FIJI/ImageJ Installation

### A. Download and installation

**If you already have FIJI, skip to step B.**

#### 1. Download FIJI:

- Website: <https://imagej.net/software/fiji/downloads>
- Choose the "Stable download" for your OS

#### 2. Installation:

- Extract the archive to a folder (e.g., `C:\Fiji` or `/Applications/Fiji.app`)
- Launch FIJI to verify installation

### B. Initial configuration

1. In FIJI: `Help` → `Update...`
  2. Click "Manage update sites"
  3. **Activate the following sites** (check boxes):
    - ☒ IBMP-CNRS
    - ☒ ImageScience
    - ☒ PTBIOP
  4. Click `Close` → `Apply changes`
  5. **Restart FIJI** after installation
- 

## STEP 2: Python and Cellpose Installation

### A. Miniconda Installation (recommended)

#### 1. Download Miniconda:

- Website: <https://repo.anaconda.com/miniconda/>
- **Windows:** `Miniconda3-latest-Windows-x86_64.exe`
- **macOS:** `Miniconda3-latest-MacOSX-x86_64.pkg`
- **Linux:** `Miniconda3-latest-Linux-x86_64.sh`

## 2. Installation:

- **Windows:** Double-click executable, check "Add to PATH"
- **macOS:** Double-click the .pkg
- **Linux:** `bash Miniconda3-latest-Linux-x86_64.sh`

## 3. Restart your terminal/command prompt

## B. Conda Configuration and Cellpose Environment Creation

### Open a terminal/command prompt and execute in order:

```
# Initial Conda configuration
conda config --set auto_activate_base false
conda config --set channel_priority flexible
conda config --set solver libmamba

# Accept Conda terms
conda info

# Clean existing environments
conda env remove -n cellpose -y
conda clean --all -y

# Create cellpose environment with Python 3.8
conda create -n cellpose python=3.8 -y

# Activate environment
conda activate cellpose
```

### Installation according to your configuration:

#### For GPU (Windows with NVIDIA):

```
pip install cellpose==3.1.1.2
pip install torch==1.13.1+cu118 torchvision==0.14.1+cu118
torchaudio==0.13.1 --index-url https://download.pytorch.org/whl/cu118 --no-
cache-dir
```

#### For CPU only (macOS/Linux or Windows without GPU):

```
pip install cellpose==3.1.1.2
pip install torch==1.13.1+cpu torchvision torchaudio --index-url
https://download.pytorch.org/whl/cpu --no-cache-dir
```

### Common dependencies:

```
pip install numpy==1.24.3 opencv-python-headless scikit-image imageio
matplotlib scipy numba
```

## C. Installation Verification

```
# Basic test
python -c "import cellpose; print('Cellpose successfully installed!')"

# GPU test (Windows with NVIDIA only)
```

```
python -c "import torch; print('CUDA available:',
torch.cuda.is_available()); print('GPU detected:',
torch.cuda.get_device_name(0) if torch.cuda.is_available() else 'None')"
```



# Cellpose GPU test

```
python -c "from cellpose import models; model = models.Cellpose(gpu=True);
print('Cellpose GPU: OK' if model.gpu else 'GPU not available')"
```

---

## STEP 3: Cellpose Plugin Configuration in FIJI

### Plugin configuration

1. In FIJI: open any image
  2. Go to Plugins → BIOP → Cellpose/Omnipose → Cellpose...
  3. Check in "to add more parameters:"
    - o  That `--use_gpu` is written (if not, add it)
    - o  That `--do_3d` is NOT present (if yes, delete it)
- 

## STEP 4: FLASH Macro Installation

### A. Macro download

1. Save the macro code in a file named `FLASH.ijm`
2. Place the file in FIJI's plugins folder:
  - o **Windows:** `C:\Fiji\plugins\`
  - o **macOS:** `/Applications/Fiji.app/plugins/`
  - o **Linux:** `/path/to/fiji/plugins/`

### B. Installation in FIJI

1. **Restart FIJI**
  2. The macro appears in Plugins → FLASH
  3. Alternative: Plugins → Macros → Run... → select `FLASH.ijm`
- 

## STEP 5: Image Preparation

### A. Required format

- **Format:** .tif or .tiff only
- **Type:** Multi-channel images (3-5 channels recommended)
- **Resolution:** Minimum 1024×1024 pixels
- **Depth:** 8-bit or 16-bit

## B. Channel organization

### Standard organization example:

1. **Channel 1:** Type IIb (MYH4) - optional
2. **Channel 2:** Laminin (membranes) - **required for segmentation**
3. **Channel 3:** Type I (MYH7) - optional
4. **Channel 4:** SDH (Oxidative) - optional
5. **Channel 5:** Type IIa (MYH2) - optional

## C. Folder preparation

1. Create a folder for analysis
  2. Copy all images to this folder
  3. Verify all files are in .tif format
  4. **Ensure channel organization is identical for all images**
- 

## STEP 6: Using FLASH

### A. Analysis launch

1. Open FIJI
2. Launch macro: `Plugins → FLASH`
3. **Fill parameters:**
  - **Total number of channels:** 4 (example)
  - **Laminin channel:** 2 (example)
  - **Fiber type channels** according to your staining
  - ☒ Check "Auto-calibrate diameter" (recommended)

### B. Advanced parameters (optional)

- **Cellpose sensitivity:** 1.0 (default)
- **Exclusion threshold:** 200 pixels
- **Adaptive threshold factor:** 0.3
- **GPU:** Check if NVIDIA card available and CUDA installed

### C. Folder selection and processing

1. Select folder containing your images
  2. Wait for processing completion (progress bar)
  3. Check results in source folder
-

## STEP 7: Results Interpretation

### A. Automatically generated files (per image)

- `ImageName_Final.tif` - Image with colored overlay
- `ImageName_Classified_Results.csv` - Classification results
- `ImageName_ROI_Set.zip` - Regions of interest
- `ImageName_cellposeMask.tif` - Segmentation mask

### B. Results analysis

- Open CSV file in Excel/LibreOffice
  - **Important columns:**
    - **ROI:** Fiber number
    - **CSA/MinFeret:** Cross-sectional area/MinFeret
    - **Classification:** Identified fiber type
    - **Fiber type statistics** at bottom of file
- 

## Troubleshooting

### Common problems and solutions

#### **CUDA Error:**

- Check graphics card compatibility
- Reinstall CUDA Toolkit 11.8
- Restart computer

#### **Cellpose not working:**

- Reinstall cellpose if necessary

#### **Images not detected:**

- Check .tif/.tiff format
- Verify channel organization
- Ensure laminin channel is correct

#### **Memory issues:**

- Close unnecessary applications
  - Use smaller image
  - Consider increasing virtual memory
-

## Additional Notes

### Performance optimization tips

- Use GPU acceleration when available for faster processing
- Process images in smaller batches if memory is limited

### Best practices

- Always backup original images before processing
- Document channel organization for each experiment
- Validate results on a subset of images before batch processing
- Keep detailed records of analysis parameters used

---

**Protocol tested on:** Windows 10/11, macOS 12+, Ubuntu 20.04+

**Version:** FLASH v3.1 – July 2025

---